# KERAMIT NP

Lega dentale Co - Cr per metallo ceramica - Tipo 4  
Co - Cr based dental alloy for metal-ceramic - Type 4  
Co - Basis Aufbrennlegierung - Typ 4  
Alliage Co - Cr pour céramo-métallique - Type 4  
Aleación Co - Cr para metalocerámica - Tipo 4

Co	Cr	W	Mo	Si	Nb
63%	24%	8%	3%	1%	1%

ISO 22674 - ISO 9693-1

CE  
0546

## Caratteristiche tecniche • Technical data • Technische daten Caractéristiques techniques • Características técnicas

Densità - Density - Dichte - Densité - Densidad	g/cm <sup>3</sup>	8,3
Intervallo di fusione - Melting range - Schmelzintervall - Intervalles de fusion - Intervalo de fusión	°C	1305 - 1370
Temperatura di colata - Casting temperature - Gießtemperatur - Température de coulée - Temperatura de colado	°C	1470
CET - CTE - WAK - CET - CET 25-500 °C - 25-600 °C	10 <sup>-6</sup> K <sup>-1</sup>	13,9 - 14,0
Limite elastico 0,2% - Yield strength 0,2% - 0,2% Dehngrenze - Limite élastique 0,2% - Limite elástico 0,2%	MPa	490
Allungamento - Elongation - Bruchdehnung - Allongement - Alargamiento	%	10
Modulo elastico - Modulus of elasticity - Elastizitätsmodul - Module d'élasticité - Módulo de elasticidad	GPa	210
Durezza Vickers - Vickers hardness - Vickershärte - Dureté Vickers - Dureza Vickers	HV10/30	285

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KERAMIT NP is Ni and Be free, it is a biocompatible non precious dental alloy suitable for the porcelain-fused to metal technique and fulfils the ISO 22674 and ISO 9693-1 standards

### INSTRUCTIONS FOR USE

- Spruing:** the wall thickness in wax has not to be thinner than 0,40 mm. Use sprues of  $\varnothing$  2,5 mm for single crowns and 3,0 mm for bridges. For prosthetic frames with more than 3 units, use a feeding wax bar of  $\varnothing$  4,0 - 4,5 mm.
- Investment:** use PVF Nobil Metal or other phosphatic investment. Pre-heating of the cylinder to 850/900°C for min. 30 minutes. Refer to the manufacturer's instructions.
- Casting:** use an individual ceramic crucible; clean the crucible after every casting to prevent from contamination with other alloys and eliminate the residuals. We recommend to pre-heat the crucible in order to shorten the casting time. If casting by induction, heat and cast when the ingots collapse and show a uniform aspect of the surface. If casting by torch (propane-oxygen) heat the metal by rotary movement of the flame and cast as soon as the metal starts vibrating. Casting temperature: 1470°C. Cool the crucible at room temperature, then free the casting from the refractory.

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**Surface preparation:** use burs in tungsten carbide specific for non precious alloys. Sand blast with aluminium dioxide 250µm at 3 bars. Clean with steam or ultrasounds.

**Oxidation:** oxidise at 960°C for 5 minutes in vacuum. Reduce the exceeding oxide by sand blasting with aluminium dioxide 250µm and further cleaning with steam or ultrasounds.

**Opaque:** we recommend, for the first layer a bonder as NM Bonding. Use a ceramic compatible with non precious alloy such as DIVA (Nobil Metal) and refer to the ceramic manufacturer's instructions.

**Finishing:** after the ceramic veneering, finish the metal parts with materials suitable for non precious alloys.

**Pre-solder:** SOLDER 1130 / KERAMIT NP SOLDER

**WARNINGS:** the metal powder is dangerous to the health. The Co-Cr based alloys can rarely cause dermatitis by contact or allergies in sensitive subjects. Use suitable aspiration systems and individual protection devices during the working phases. Before applying the prosthetical medical device, verify if other metal implants are already existing in the oral cavity of the patient to prevent galvanic phenomena.

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